EN Mediblink Infrared Forehead Thermometer Age Accuracy™ M330

INSTRUCTIONS FOR USE

Please read the instructions for use carefully before using the product

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Pažljivo pročitajte upute prije upotrebe uređaja

TABLE OF CONTENTS

Introduction	5
Warning	6
Product Description	7
LCD Display Introduction	8
Basic Functions	9
Temperature Taking Hints	10
Illustration For Use	11
Memory Mode	14
Care And Cleaning	15
Battery Replacement	16
Technical Specifications	17
Troubleshooting	18
Calibration	20
Symbol Explanation	21
Electromagnetic Compatibility Information	22
Warranty	28

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Introduction

Please read all instructions carefully and thoroughly before using this product.

Infrared forehead thermometer Age Accuracy™ M330 (DET-3017) is specifically designed for safe use on the forehead. The Infrared Forehead Thermometer is a device capable of measuring people's body temperature by detecting the intensity of infrared light emitted from the forehead. It converts the measured heat into a temperature reading displayed on the LCD. When properly used, it will quickly assess your temperature in an accurate manner.

Intended use: The infrared forehead thermometer is intended for the intermittent measurement of human body temperature from the skin surface of forehead. The device can be reused by people of all ages for home use and clinical use.

Indications for Use: The infrared forehead thermometer is used to measure body temperature by measuring forehead.

This appliance conforms to the following standards:

ISO 80601-2-56 Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement,

IEC 60601-1-11 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of IEC 60601-1-2(EMC),



IEC 60601-1(Safety) standards. And the manufacturer is ISO 13485 certified.

Warning


1. There is no gender or age limitation for using infrared forehead thermometer.
2. Do not touch the temperature probe with hands.
3. Use of this Forehead thermometer is not intended as a substitute for consultation with your physician.
4. Do not allow children to take their temperatures unsupervised, some parts are small enough to be swallowed.
5. Never immerse this device in water or other liquids.
6. Do not modify this equipment without authorization of manufacturer.
7. Do not expose the thermometer to temperature extremes (below -25°C/-13°F or over 55°C/131°F) nor excessive humidity (>95%RH).
8. Keep the battery away from children.
9. Remove battery from the device when not in operation for a long time.
10. Do not put the thermometer in direct sunlight or with cotton wool, otherwise the accuracy will be affected.
11. Portable and mobile RF communications can affect the device. The device needs special precaution regarding EMC according to the EMC information provided in the accompanying documents.
12. ME equipment should not be cleaned while in use.
13. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.
14. The probe of the ME equipment shall not be serviced or maintained while in use with a patient.

Product Description

1. Probe
2. Test Button (for measuring)
3. Sound Switch
4. Setting Button
5. Battery Cover
6. Display

LCD Display Introduction

Basic Functions

Forehead Mode	The thermometer has been designed for practical use. It's not meant to replace a visit to the doctor. Please also remember to compare the measurement result to your regular body temperature. → Please see the Illustration For Use section to learn how to measure the body temperature.
Object Mode	The object mode shows the actual, unadjusted surface temperature, which is different from the body temperature. It can help you to monitor if the object temperature is suitable for the baby or patient, for example the baby's milk. Measuring range of object mode: 0°C~100°C (32°F~212°F) Laboratory accuracy of object mode: ±4% or ±2°C (4°F) whichever is greater. → Please see the Illustration For Use section to learn how to measure the object temperature.
Memory Mode	There are 30 sets memories for forehead and object measurements. Each memory also records the measurement mode and age range icon.
°C / °F Switch	Please see the Illustration For Use to learn how to change between Celsius and Fahrenheit.
Beep Hint	When the testing is done, it will sound one short beep.
Sound on/off	The thermometer can turn on/off <u>sound</u>. → Please see the Illustration For Use.
Age range	Please see Illustration For Use section to learn how to select age range <u>among</u> 0–3 months, 3–36 months and 36 months–adult.

Temperature Taking Hints

To ensure that the reading always reflects the body temperature accurately, you need to take account of the following factors which may affect an accurate reading.

1. It is important to know each individual's normal temperature when they are well. This is the only way to accurately diagnose a fever. To determine normal temperature, take multiple readings when healthy. Re-measure with a standard digital thermometer for confirmation.
2. Users must be inside for 30 minutes before taking a measurement. Note: Users and the thermometer should be in the same ambient temperature for at least 10 minutes before taking a reading.
3. Users should not drink, eat, or be physically active such as bathing, showering, shampooing and hair drying before/while taking the measurement. Remove hat and hair and wait 10 minutes before taking a reading.
4. Oils or cosmetics on the forehead may give a lower temperature reading than the actual one. Remove dirt from the forehead before taking a measurement. Wait at least 10 minutes after washing the forehead area before taking a reading.
5. Holding a hand on the forehead for any length of time will affect the temperature reading.
6. Do not take temperature over scar tissue, open sores or abrasions.
7. Do not use the thermometer on a perspiring or sweating forehead, as this may affect the reading.
8. Don't take a measurement while or immediately after nursing a baby.
9. Do not take temperatures with this thermometer near places that are very hot, such as fireplaces and stoves.
10. The probe window of the thermometer is the most delicate part of the device. Do not touch the probe window. The accuracy of the reading may be affected if the probe window is damaged or dirty.

11. If the thermometer is stored in a significantly different environment than testing location, place it in the testing location for approximately 30 minutes prior to use.
12. It is not intended for use in the oxygen rich environment and presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
13. The patient is an intended operator.

Illustration For Use

➤ To measure forehead temperature:

1. Press the *Test Button*, The display is activated to show all segments. After self-checking Figure 1 appears on the display screen with beep, so you can start a new measurement.

Note: You can press *Setting Button* to change age range among 0–3 months, 3–36 months and 36 months–adult.

2. Aim the thermometer at the center of the forehead with a distance less than 5cm (See figure 2) and then press the *Test Button*.

Note: Do not remove the thermometer from the forehead before hearing beep.

3. Read the temperature on the display.
4. Device will automatically shut off if left idle for 30 seconds.


Figure 1

Figure 2

➤ How to turn on or off sound:

You can flip the *Sound Switch* to turn on or off sound.

➤ How to change the forehead mode and object mode:

You can press and hold *Setting Button* to switch the mode between Forehead mode and Object mode.

➤ How to set measurement unit between °C and °F:

With the thermometer off, press and hold *Setting Button* to enter into setting mode. Press *Test Button* to select the unit you want. After the unit is set, press *Setting Button* to exit the setting mode (See Figure 3).


Figure 3

➤ How to select the age range:

With the thermometer on, press the *Setting Button* to select the age range among 0–3 months, 3–36 months and 36 months–adult (See Figure 4).


Figure 4

(Default)

Age range	GREEN	YELLOW	RED
 0–3 Months	34.0–37.4°C (93.2–99.4°F)	×	37.5–43.0°C (99.5–109.4°F)
 3–36 Months	34.0–37.6°C (93.2–99.6°F)	37.7–38.5°C (99.7–101.3°F)	38.6–43.0°C (101.4–109.4°F)
 36 Months–adult	34.0–37.7°C (93.2–99.9°F)	37.8–39.4°C (100.0–103.0°F)	39.5–43.0°C (103.1–109.4°F)

➤ To measure object temperature:

1. Press the *Test Button* to turn on the thermometer, you can take the object temperature after going into measurement mode (See Figure 5).
2. Aim the thermometer at the center of the object you want to measure with a distance less than 5cm.
3. Press the *Test Button* and then read the temperature on the display.
4. Device will automatically shut off if left idle for 30 seconds.


Figure 5

➤ After measurement:

1. Power off: Device will automatically shut off if left idle for more than 30 seconds to extend battery life.
2. Clean the probe after each use to ensure an accurate reading and avoid cross contamination.
(See the section of Care and Cleaning for details)

> BackLight:

In Forehead mode:

The display will be lighted for 3 seconds when a measurement is completed.

In Object mode:

The display will only be lighted WHITE for 3 seconds when the unit is ready for measurement and a measurement is completed.

Memory Mode

1. With the thermometer off, press the *Setting Button*. The letter M will appear in the upper left corner of the display. (See Figure 6)
2. The thermometer will automatically memorize the last 30 temperature readings. Each memory also records the measurement mode icon. Each time the *Setting Button* is pressed, the screen displays past readings that correspond with a number 1–30. The number 1 reflects the most recent reading, while the number 30 reveals the oldest reading stored in memory. (See Figure 7)
3. In the memory mode,  mark or  mark or range mark will not change. The user can press the *Test Button* to take new measurements.


Figure 6


Figure 7

Care And Cleaning

1. The probe window must be kept clean, dry, and undamaged at all times to ensure accurate readings. The accuracy of temperature readings can be affected by damage to the probe window, or the presence of dirt, fingerprints, dust and other soiling compounds on the probe window. Degraded sensors can degrade performance or cause other problems.
2. The probe window is the most delicate part of the thermometer. Use a soft cloth slightly moistened with a 75% isopropyl alcohol solution to wipe probe window and the thermometer. Do not use abrasive cleaners. After cleaning, allow at least 10 minutes drying time before taking temperatures.
Note: Do not use any chemical other than isopropyl alcohol to clean the probe window.
3. Use a soft, dry cloth to clean the thermometer display and exterior.
4. Do not put the thermometer into water directly.
5. Store the thermometer in a dry location, free from dust and contamination and away from direct sunlight.
6. Put the thermometer back to the original packaging after using.

Battery Replacement

1. Replace battery when “” flashes in the upper right corner of LCD display. (See Figure 8)
2. Slide battery cover down as shown in Figure 9.
3. Remove battery and install 2 new AAA alkaline batteries as shown in Figure 10.
4. Slide battery cover back on.


Figure 8


Figure 9


Figure 10


Technical Specifications

Model / Reference number	DET-3017 / M330
Measuring range	Forehead mode: 34.0~43.0°C (93.2~109.4°F)
Measuring site	Forehead (Forehead Mode)
Reference body site	Oral (This thermometer converts the forehead temperature to display its “oral equivalent.”)
Operation mode	Forehead mode (Adjust mode)
Laboratory accuracy	Forehead mode: ±0.2°C (0.4°F) during 35.5~42.0°C (95.9~107.6°F) at 15~35°C (59.0~95.0°F) operating temperature range ±0.3°C (0.5°F) for other measuring and operating temperature range
Display resolution	0.1°C or 0.1°F
Measure time	Approximately 1 second
Operating temperature range	5~40°C (41~104°F), 15~85%RH, non-condensing Atmospheric Pressure : 70~106kPa
Storage and transport temperature range	-25~55°C (-13~131°F), 15~95%RH, non-condensing Atmospheric Pressure: 70~106kPa
Clinical accuracy	0–1 year: Clinical bias: -0.19°C (-0.34°F); Clinical repeatability: 0.06°C (0.11°F); Limits of agreement: 0.68°C (1.22°F) 1–5 years: Clinical bias: -0.25°C (-0.45°F); Clinical repeatability: 0.07°C (0.13°F); Limits of agreement: 0.73°C (1.31°F) over 5 years: Clinical bias: -0.25°C (-0.45°F); Clinical repeatability: 0.17°C (0.31°F); Limits of agreement: 0.60°C (1.08°F)

Shock	Withstands drop of 3 feet
Dimension	155*39.2*49.4mm
Weight	Approx. 97 grams (with batteries)
Battery	DC3V (2×AAA battery)
Battery life	Approx. 6000 readings
Applied part	Probe
Expected service life/ Shelf life	Three years
Ingress protecting rating	IP22
Accessories	Battery
Contraindications	None

Troubleshooting

Error message	Problem	Solution
Er 1	The thermometer is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
Er 2	The ambient temperature is not within the range between 5°C and 40°C (41°F~104°F).	Place the in a room for at least 30 minutes at room temperature between 5°C and 40°C (41°F~104°F).

Error message	Problem	Solution
	In Forehead mode: Temperature taken is higher than 43.0°C (109.4°F). In Object mode: Temperature taken is higher than 100°C (212°F).	Read Temperature Taking Hints Thoroughly, then take a new temperature measurement.
	In Forehead mode: Temperature taken is lower than 34.0°C (93.2°F). In Object mode: Temperature taken is lower than 0°C (32°F).	Read Temperature Taking Hints thoroughly, then make sure the lens filters are clean, then take a new temperature measurement.
	The thermometer works properly.	Use the thermometer normally
	When voltage sign flashes, it indicates that the power is low, but you can continue to measure.	The thermometer will take a proper measurement but batteries must be replaced soon.
	When voltage sign flashes with nothing else, the thermometer can not work due to low battery.	Replace with two new alkaline batteries size AAA.

Calibration

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the instructions for use, periodic readjustment is not required. However, We recommend checking calibration every year or whenever clinical accuracy of the thermometer is in question. Please send the complete device to the dealer or manufacturer.

The above recommendations do not supersede the legal requirements. The user must always comply with legal requirements for the control of the measurement, functionality, and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

A clinical summary and procedures for checking calibration are available upon request (Turn on the thermometer and press the setting button long time until entering into calibrate mode, software version will be displayed).

Symbol Explanation

	Manufacturer		Storage and Transportation Humidity limitation: 15%~95%RH
	Authorised Representative in EU		Electrical and electronic equipment marking. Reduce electronic and electrical waste as unsorted waste and collect it separately.
	Date of manufacture	IP22	The first num.2: Protected against solid foreign objects of 12,5 mm $\varnothing$ and greater. The second num.2: Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
	Product reference number		
	Batch number*		
	Unique Device Identifier		CE Marking of Conformity, Notified body ID number
	Storage and Transportation Temperature Limit: -25°C~55°C (-13°F~131°F)		Medical Device
	Type BF applied part		Caution
	Before use, read the instructions. Electronic instructions for use: http://www.mediblink.com/f/m330.pdf		
	Atmospheric pressure limitation		

*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

Electromagnetic Compatibility Information

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1:

Guidance and manufacturer's declaration – electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
Conducted emission CISPR 11	Not applicable	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radiated emission CISPR 11	Group 1 Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	
		The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Table 2:

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	±2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	N/A	N/A
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode line-line	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T (100% dip in U_T) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T (100% dip in U_T) for 1 cycle at 0° 70% U_T (30% dip in U_T) for 25/30 cycles at 0° 0% U_T (100% dip in U_T) for 250/300 cycle at 0°	N/A	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity magnetic fields	IEC 61000-4-39	See Table 3	See Table 3
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Table 3:

Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields		
Test frequency	Modulation	IMMUNITY TEST LEVEL (A/m)
30 kHz ^a	CW	8
134,2 kHz	Pulse modulation ^b 2,1 kHz	65 ^c
13,56 MHz	Pulse modulation ^b 50 kHz	7,5 ^c

^aThis test is applicable only to ENVIRONMENT. ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE

^bThe carrier shall be modulated using a 50% duty cycle square wave signal.

^cr.m.s., before modulation is applied.

Table 4:

Guidance and manufacturer's declaration – electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands ^a	N/A
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
^a The ISM(industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHzto6,795 MHz;13,553 MHz to 13,567 MHz;26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz,3,5 MHz to 4,0 MHz,5,3 MHz to 5,4 MHz,7 MHz to 7,3 MHz,10,1 MHz to 10,15 MHz,14 MHz to 14,2 MHz,18,07 MHz to 18,17 MHz,21,0MHz to 21,4MHz,24,89 MHz to 24,99 MHz,28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.		
^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.		
^c Field strengths from fixed transmitters, such as base stations for radio(cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.		
^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Table 5:

Recommended separation distances between RF wireless communications equipment				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.				
Frequency (MHz)	Maximum Power (W)	Distance (m)	IEC 60601 Test Level	Compliance Level
385	1.8	0.3	27	27
450	2	0.3	28	28
710	0.2	0.3	9	9
745				
780				
810	2	0.3	28	28
870				
930				
1720	2	0.3	28	28
1845				
1970				
2450	2	0.3	28	28
5240	0.2	0.3	9	9
5500				
5785				
Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

 WARNING

This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.

- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Portable and mobile RF communications can affect the device. The device needs special pre-cautions regarding EMC according to the EMC information provided in the accompany documents.
- Do not use the devices in the MR environment.
- The Operator should not use the system and should inform the customer service, if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES.
- PRECAUTION: The performance of the device may be degraded should one or more of the following occur:
 - Operation outside the manufacturer's stated temperature and humidity range.
 - Storage outside the manufacturer's stated temperature and humidity range.
 - Mechanical shock (for example, drop test) or degraded sensor.
 - Patient temperature is below ambient temperature.

Warranty

Product: Mediblink Infrared Forehead Thermometer Age Accuracy™ M330

Manufactured for (Importer for EU & Distributor): Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

Seller's name, address, signature and stamp*:

Date of extradition / sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

WARRANTY TERMS

Dear customers!

The warranty period is **3 years** and starts from the day of product purchase. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product cannot be repaired or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product cannot be repaired in 45 days, the product will be replaced with a new one. In case that the product cannot be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Manufacturer:
JOYTECH Healthcare Co., Ltd.
No.365, Wuzhou Road,
311100 Hangzhou, Zhejiang Province,
PEOPLE'S REPUBLIC OF CHINA
info@sejoy.com

Manufactured for
(Importer for EU & Distributor):
Mediblink d.o.o.
Gubčeva cesta 19
8210 Trebnje
Slovenia
info@mediblink.com
www.mediblink.com

Shanghai International Holding Corp.
GmbH (Europe)
Eiffestrasse 80,
20537 Hamburg,
Germany
shholding@hotmail.com

Instructions for use version: V.01
Issue date: 15. 7. 2024
Date of last change: 25. 3. 2025

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